


## Accessories

---

### Baskets

Use of a basket is recommended. Baskets prevent objects touching the sides of the tank and allow items to be placed in or removed from the bath without immersing your hands in the chemical solution.

<b>Bath</b>	<b>Suitable basket</b>	<b>Dimensions l/w/d mm</b>	<b>Overall depth to tank rim mm</b>
XB2	AB2	102/89/64	93
XB3	AB3	195/120/64	93
XB6	AB6	237/102/102	130
XB14	AB14	263/236/102	130
XB22	AB22	445/252/102	130



### Lids

Lids minimise operating noise and reduce evaporation

<b>Bath</b>	<b>Suitable lid</b>
XB2	AL2
XB3	AL3
XB6, MXB6	AL6
XB14, MXB14	AL14
XB22, MXB22	AL22